

Optomotry

DESCRIPTION


Optomotry is available for the measurement of visual acuity in mice. Animals are placed in the apparatus with a rotating black and white gradient. The experimenter determines whether the animal appears to be following the rotation of the gradient. The Animal Phenotyping Core houses one of these units.

FEATURES & SPECIFICATIONS

Four monitors facing inward towards the apparatus display virtual stimuli

Software allows continuous re-centering of the virtual stimuli based on the animal's position, programming of experiments, system calibration, and data logging

PHOTOS



Our Mission:
Dedicated to the work of
health and healing