

## Visual Cliff

### DESCRIPTION

The Visual Cliff test is available for the measurement of visual acuity. Animals are placed in a stable apparatus with a clear Plexiglas bottom. The apparatus is placed over the ledge of a table, and the amount of time the animal spends over the 'visual cliff' is recorded. The Animal Phenotyping Core houses two of these units for use on mice, rats, and other small animal models.

### FEATURES & SPECIFICATIONS


Stable wooden frame with clear Plexiglas bottom

Coupled with Any-Maze behavioral tracking software for objective data collection and analysis

Any-Maze behavioral tracking software

- Version 5.2

### PHOTOS



Our Mission:  
Dedicated to the work of  
health and healing